

# Hydrofera Blue CLASSIC

Polyvinyl Alcohol (PVA) Foam Dressings

# Hydrofera Blue READY

Polyurethane (PU) Foam Dressings

Antibacterial foam dressings that combine methylene blue and gentian violet to manage bioburden.

## Wound characteristics

- Deep, tunneling or shallow
- With or without devitalized tissue/slough & debris
- Moderate to heavy exudate
- Suspected biofilm
- Infected or at risk of infection

## Clinical benefits

- Continuous wicking, highly absorptive
- Safe, autolytic debridement
- Disrupts biofilm
- Manages bioburden
- Safe, Non-cytotoxic
- Does not inhibit growth factors
- High tensile strength



## Wound characteristics

- Shallow wounds
- Moderate to heavy exudate


## Clinical benefits

- Up to 7-day wear time
- Non-cytotoxic
- Does not inhibit growth factors



## Powerful Wicking

- Wicks bacteria laden exudate into the foam and away from the wound surface with natural negative pressure through capillary flow.



## Debridement

- Example of devitalized tissue removed at dressing change.



## Edge Effect

- Multiple authors have noted resolution of epibole or flattening of wound edges with use of Hydrofera Blue CLASSIC dressings.<sup>1-4</sup>



Granulating wound bed with epithelial cells on the edge

## Disrupts Biofilm

- A recent meta-analysis of nine studies involving 185 chronic wounds found that 78.2% of those wounds contained bacterial biofilm.<sup>5</sup>
- Biofilms in wounds are known to form on, and associate with slough and devitalized tissue.<sup>6,7</sup>
- Hydrofera Blue CLASSIC dressing manages exuding wounds and aids in the removal of devitalized tissue from the wound bed<sup>8-11</sup> As this action occurs, Hydrofera Blue CLASSIC dressing helps disrupt biofilm that is associated with slough and devitalized tissue.



## Film Backing

- Waterproof and bacteria-resistant semi-permeable backing


## Broad-spectrum antibacterial protection

## Absorption


- Absorbs and retains bacteria-laden exudate



Deep with Heavy Exudate



Tunneling



Removal of Devitalized Tissue



Shallow with Moderate Exudate



Skin Tear Before / After



With Compression / Total Contact Casting

Safe, Effective, Protection from Start to Finish

All Hydrofera Blue dressings are effective against the broad-spectrum of microorganisms commonly found in wounds.

Bacteria			Yeast
<ul style="list-style-type: none"> <li>Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)</li> <li><i>Staphylococcus aureus</i></li> <li><i>Staphylococcus epidermidis</i> (coagulase negative)</li> <li><i>Pseudomonas aeruginosa</i></li> </ul>	<ul style="list-style-type: none"> <li><i>Pseudomonas fluorescens</i></li> <li><i>Escherichia coli</i></li> <li><i>Enterococcus faecalis</i></li> <li>Vancomycin-resistant <i>Enterococcus</i> (VRE)</li> </ul>	<ul style="list-style-type: none"> <li><i>Streptococcus pyogenes</i></li> <li><i>Klebsiella pneumoniae</i></li> <li><i>Proteus mirabilis</i></li> <li><i>Proteus vulgaris</i></li> <li><i>Enterobacter aerogenes</i></li> </ul>	<ul style="list-style-type: none"> <li><i>Bacillus subtilis</i></li> <li><i>Yersinia enterocolitica</i></li> <li><i>Serratia marcescens</i></li> <li><i>Streptococcus pneumoniae*</i></li> <li><i>Enterobacter cloacae*</i></li> </ul>
<ul style="list-style-type: none"> <li><i>Candida albicans</i></li> <li><i>Candida krusei</i></li> <li><i>Candida glabrata</i></li> </ul>			

\* Tested on Hydrofera Blue Ready dressing only

## ORDERING INFORMATION - Hydrofera Blue Antibacterial Foam Dressings

	Product No	Product Size	Qty/Box
<p>No hydration required</p> <p>Hydrate before use</p>	<b>Hydrofera Blue CLASSIC</b>		
	HB2214	5 cm x 5 cm (2" x 2")	10
	HB4414	10.2 cm x 10.2 cm (4" x 4")	10
	HB6614	15.2 cm x 15.2 cm (6" x 6")	10
	HBT0906	9 mm Tunneling Dressing (1.2 g)	10 pair
	HBRF2650	6.4 cm (2.5") diameter Ostomy Dressing	10
<p>No hydration required</p>	<b>Hydrofera Blue CLASSIC Heavy Drainage</b>		
	HBHD4450	10.2 cm x 10.2 cm x 1.27 cm (4" x 4" x 0.5")	5
	HBHD6675	15.2 cm x 15.2 cm x 1.9 cm (6" x 6" x 0.75")	5
<p>No hydration required</p>	<b>Hydrofera Blue READY</b>		
	HBRS2520	6.4 cm x 6.4 cm (2.5" x 2.5")	10
	HBRS4520	10.2 cm x 12.7 cm (4" x 5")	10
	HBRS8820	20.3 cm x 20.3 cm (8" x 8")	10

Indications for use: Pressure Injuries/Ulcers, Diabetic Ulcers, Venous Stasis Ulcers, Arterial Ulcers, Superficial Burns, Donor Sites, Skin Grafts, Post Surgical Incisions, Traumatic Wounds, Skin Tears, Abrasions, Lacerations.

### REFERENCES

- Edwards K. New twist on an old favorite: gentian violet and methylene blue antibacterial foam dressings. *Adv Wound Care*(New Rochelle). 2016; Jan 1;5(1):11-18.
- Swan H, Trovela VJ. Case study review: use of an absorbent bacteriostatic dressing for multiple indications. Poster presented at Clinical Symposium on Advances in Skin and Wound Care; September 9-11, 2011; Washington, DC.
- Weir D, Blakely M. Case review of the clinical use of an antimicrobial PVA foam dressing. Poster presented at Symposium on Advances in Skin and Wound Care; April 18-21, 2012; Atlanta, GA.
- Conwell P, Mikulski L, Tramontozzi M. A comparison of two antimicrobial PVA foam dressings: a randomized prospective trial comparing PVA foam with two organic pigments to a silver based wound dressing. Poster presented at Symposium on Advanced Wound Care, May 2-5, 2004; Lake Buena Vista, Fla.
- Malone M, Bjarnsholt T, McBain AJ, et al. The prevalence of biofilms in chronic wounds – a systematic review and meta-analysis of published data. *J Wound Care*. 2017; Jan 2;26(1):20-25.
- Percival SL, Suleman L. Slough and biofilm: removal of barriers to wound healing by desloughing. *J Wound Care*. 2015; Nov;24(11):498-510.
- Nakagami G, Schultz G, Gibson DJ, et al. Biofilm detection by wound blotting can predict slough development in pressure ulcers: a prospective observational study. *Wound Rep and Reg*. 2017; 25:131-138.
- Applewhite AJ, Attar P, Liden B, Stevenson Q. Gentian violet and methylene blue polyvinyl alcohol foam antibacterial dressing as a viable form of autolytic debridement in the wound bed. *Surg Technol Int*. 2015 May; 26:65-70.
- Hill R. Optimizing the wound bed by removing devitalized tissue and using methylene blue and gentian violet antibacterial foam dressings: a case series. Poster presented at Wounds Canada; May 12-14, 2017; Kamloops, BC.
- Prest D. Managing challenging chronic wounds in the community setting using an antibacterial PVA foam dressing containing methylene blue and gentian violet. Poster presented at CAWC; October 29 – November 1, 2015; Toronto, ON.
- Woo KY, Heil J. A prospective evaluation of methylene blue and gentian violet dressing for management of chronic wounds with local infection. *Int Wound J*. 2017; doi: 10.1111/iwj.12753

For product questions and samples please call Customer Service at 1-860-337-7730.

Hydrofera Dressings manufactured by:  
Hydrofera, LLC, 340 Progress Drive  
Manchester, CT 06042, USA  
[www.hydrofera.com](http://www.hydrofera.com)



Hydrofera Blue and the Hydrofera Blue logo are trademarks of Hydrofera, LLC. ©2018  
Hydrofera Blue Dressings are represented by Appulse.

**Hydrofera**®

# CHANGE the OUTCOME



with **Hydrofera Blue**

## Hydrofera Blue

- Non-Cytotoxic
- Natural Negative Pressure
- Natural Autolytic Debridement
- Antibacterial
- Flattens Rolled Wound Edges
- Disrupts Biofilm

Change to blue at [hydrofera.com](http://hydrofera.com)

Hydrofera**BLUE**