

Introducing CollarForm[®] Arm Sling:



SPLINTS & IMMOBILISATION



CollarForm[™]

 $CollarForm^{M}$ is a strong, stretch resistant sling made of a 5cm wide strip of dense foam covered with a pure cotton stockinette.

 $CollarForm^{M}$ is ideally designed to immobilise the clavicle.

 $CollarForm^{TM}$ is 100% Australian Made.

ode	Width	Length	Packaging
M000	5cm	10m	Economy Roll
M001	5cm	2 x 6m	Dispenser Pack
MUU1	SCM	2 x 6m	Dispenser Pack