

Botox Vacc-Safe

Cosmetic Refrigerator

BT42

LED Display

Microprocessor control, digital display and temperature adjustment with an increment of 0.1 C.

Multiple Alarm Functions

Capable of alerting failures due to high and low temperature, sensor error, power failure, low battery, door ajar and remote alarm.

Ergonomic Design

810mm cabinet height makes it easy to place under test bench, saving storage space.

Lockable Design

Safety lock to prevent unauthorized access.

USB Data Recording

Easy use USB for recording temperature data



Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



Audible and Visual Alarms

Audible buzzer and visible flashing light



Interior
LED



Superior Cabinet Temperature Uniformity, Auto Defrost

0.1°C Temperature accuracy

Energy Efficient Hydrocarbon (HC) Refrigeration System

Optimised refrigeration system, utilising high-efficiency HC compressor significantly reduces energy consumption.

Specifications

Technical Data	
Cabinet Type	Upright, glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	R600a
Noise (dB (A))	35
Performance	
Temp Range (C)	2~6

Construction		
Capacity (L)		42
Net/Gross Weight	(kg)	25/29
Interior Dimension (W*D*H) (mm)		400*345*380
Exterior Dimension	(mm)	480*450*550
(W*D*H) Packing	(mm)	540*640*710
Electrical Data		
Power Supply (V/Hz)		220-240/50 (60)
Power (W)		65
Electrical Current (A)		0.75

Alarms	
High/Low Temp	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

EBOS Code:

49520106