




## COVER-DRI

Cover-Dri is the Attends range of surface protection underpads, specially designed to manage light to moderate urinary incontinence and for clinical procedures where non-sterile containment is required.



- 1 The **backsheet material** helps to keep the product in place.
- 2 Discreetly coloured **waterproof backsheet** helps to protect surfaces.
- 3 **Hexagonal quilt design** with individual cushions containing fluff and/or super absorbents to quickly absorb and lock away fluid.
- 4 **Drainage channels** rapidly distribute the fluid without causing it to pool in one place.
- 5 **All-round leakage barrier** helps to prevent leakage on to surfaces.
- 6 **Super absorbent core** provides leakage protection, skin dryness and odour protection.
- 7 Attends Cover-Dri Plus 80 x 170cm features **tuckable wings**, which can be used on beds and other surfaces to securely keep the pad in place.



Skin tolerability is dermatologically approved by SGS proderm Institute for Applied Dermatological Research.

For more information visit [www.attends.co.nz](http://www.attends.co.nz)

**EBOS** HEALTHCARE [www.eboshealthcare.co.nz](http://www.eboshealthcare.co.nz)




**Super Soft**  
Feels soft against  
your skin.



**Odour Protection**  
Keeps you safe from  
unpleasant smells.



**Waterproof Backsheet**  
Helps to keep  
surfaces dry.

<b>Attends COVER-DRI</b> 	<b>Dimensions (cm)</b>	<b>Qty Pack Case</b>	<b>GTIN (7332152) Pack Case</b>	<b>EBOS Codes</b>
<b>Plus</b>	40 x 60	50 6 x 50	203903 203910	10003891
	60 x 60	50 4 x 50	203927 203934	10003892
	60 x 90	50 4 x 50	203965 203972	10003893
	80 x 90	30 4 x 30	204009 204016	10003894
<b>Plus (with Tuckable Wings)</b>	80 x 170	30 4 x 30	204023 204030	10003895
<b>Super</b>	60 x 90	50 4 x 50	203989 203996	10003897

### Materials

#### Absorbent Core Material

Plus/Super - Super absorbents and fluff

#### Back Sheet Material

Polypropylene film

#### Top Sheet Material

Non-woven

#### Latex free



**Attends Healthcare AB**  
Järnvägsgatan 4  
SE-57833 Aneby  
Sweden



At Attends, we care about each other,  
our community and the planet.



Products produced by using  
renewable energy and  
supporting verified emission  
reduction projects.



These products have  
achieved the certification  
of the recognised Nordic  
Swan Ecolabel.



Plastic bags made  
with min. 30%  
recycled material.



The cellulose within these  
products originates from  
sustainably managed forests  
and controlled sources.



With this label these  
products are certified  
as being harmless for  
human health.