



Comfortable
for Patient



Fast
Absorbency



Gentle
on Skin

Mediam 135 Low Adherent Dressing - Lite consists of a low profile absorbent cotton inner pad encased between two layers of high density polyethylene film.

LOW-ADHERENT

Low adherence to the wound, therefore minimal trauma and pain upon removal

ABSORBENT INNER PAD

Low profile inner pad, wicks excess exudate away from wound bed

CUT TO SIZE

Maintains integrity when cut to size

Suitable for:

Light exuding and small wounds


First aid applications

Minor cuts and grazes

EBOS Code	Reorder	Size	Sterility	Pcs/Unit
28844003	NAL001	7.5 x 5cm	Sterile	100
	NAL002	7.5cm x 7.5cm	Sterile	100
28844004	NAL003	7.5cm x 10cm	Sterile	100
	NAL004	7.5cm x 20cm	Sterile	100



FOR MORE INFORMATION CONTACT YOUR ACCOUNT MANAGER
ON **FREEPHONE 0800 105 501** OR VISIT www.eboshealthcare.co.nz



ESSENTIAL TO HEALTH