



Vernacare

VORTEX AIR

The new platform for
macerator technology

Smarter diagnostics

Sustainable design

Superior performance



ESSENTIAL TO HEALTH



Technical Specifications

The Vortex Air utilises the latest in motor control and diagnostics tools to minimise machine downtime whilst our SmartFlow™ maceration system ensures maximum performance. The smooth design and antimicrobial surfaces also work to reduce the risks of HCAI's alongside our sensor controlled Touchless Technology operation.

Cycle Time (Standard)	120 seconds
Cycle Time (Heavy Duty)	140 seconds
Cycle Time (Eco Mode)	95 seconds
Weight	76kg
Capacity	6 small, 4 medium or 2 large disposable pulp products
Closing Mechanism	Automatic (via foot sensor)
Opening Mechanism	Automatic (via foot sensor)
Outside Dimensions	427 x 628 x 1027mm
Drain Outlet	2" P Trap (included)
Water Inlet	1/2" BSP (male)
Water Usage (Standard)	23L
Water Usage (Heavy Duty)	27L
Water Usage (Eco Mode)	16L
Required Flowrate	12L/min (minimum)
Sound Level	62 dBA (average)

EBOS Code: 11000410

FOR MORE INFORMATION CONTACT YOUR ACCOUNT
MANAGER ON **FREephone 0800 105 501** OR VISIT
www.eboshealthcare.co.nz

 **EBOS**
HEALTHCARE

ESSENTIAL TO HEALTH