



NEW ENDOFORM® PRESENTATIONS NOW AVAILABLE

The **Endoform®** portfolio has expanded to provide solutions for wounds of varying exudate levels and sizes

Endoform® 1" Disc for wounds 1 inch in size or less

- Simplifies application with reduced trimming and wastage
- Single use presentation which is intended to ensure a sterile piece is applied each time
- Optimal shape for small, circular wounds e.g. diabetic foot ulcers



 Available in both Natural and Antimicrobial presentations


| Endoform® Natural - 1" Disc | | |
|-----------------------------|--------------|--------------|
| Stock no. | Product Size | Quantity/Box |
| 529315 | 1" Disc | 10 |

| Endoform® Antimicrobial - 1" Disc | | |
|-----------------------------------|--------------|--------------|
| Stock no. | Product Size | Quantity/Box |
| 629315 | 1" Disc | 10 |

Endoform® High Flow: with larger apertures for wounds with a variable level of exudate

- Can be used in combination with Hydrofera™ Blue Classic
- Ideal for use with Negative Pressure Wounds Therapy (NPWT)



 Available in both Natural and Antimicrobial presentations



| Endoform® Natural - High Flow | | |
|-------------------------------|----------------------------|--------------|
| Stock no. | Product Size | Quantity/Box |
| 529302 | 2x2" (5cm x 5cm) | 10 |
| 529304 | 4x5" (10cm x 12.7cm) | 10 |
| 529316 | 6x8.5" (15.24cm x 21.59cm) | 10 |

| Endoform® Antimicrobial - High Flow | | |
|-------------------------------------|----------------------|--------------|
| Stock no. | Product Size | Quantity/Box |
| 629302 | 2x2" (5cm x 5cm) | 10 |
| 629304 | 4x5" (10cm x 12.7cm) | 10 |

Endoform® is unique extracellular matrix (ECM) to support all phases of wound healing

It provides a natural, porous ECM scaffold for rapid cell infiltration and 148 secondary molecules important for healing.

ORDERING INFORMATION FOR OTHER ENDOFORM[®] PRESENTATIONS



Endoform[®] Natural Dermal Template - Fenestrated

| Stock no. | Product Size | Quantity/Box |
|-----------|----------------------------------|--------------|
| 529312 | 2x2" (5cm x 5cm) fenestrated | 10 |
| 529314 | 4x5" (10cm x 12.7cm) fenestrated | 10 |

Endoform[®] Antimicrobial Dermal Template - Fenestrated

| Stock no. | Product Size | Quantity/Box |
|-----------|----------------------------------|--------------|
| 629312 | 2x2" (5cm x 5cm) fenestrated | 10 |
| 629314 | 4x5" (10cm x 12.7cm) fenestrated | 10 |

RX Only. Prior to use, be sure to read the entire Instructions For Use package insert supplied with the product.

Endoform[®] is a registered trademark of Aroa Biosurgery Limited.
Hydrofera[™] Blue is a registered trademark of Hydroferra LLC.



Manufactured for:
AROA BIOSURGERY INC

2 Kingsford Smith Place,
 Auckland 2022, New Zealand
 1-877-627-6224
www.aroabio.com