



 **smith&nephew**  
**CUTIPLAST<sup>®</sup>**  
**STERIL**  
Conformable Fabric Island  
Dressing

A dressing for minor injuries

# CUTIPLAST<sup>®</sup> STERIL Conformable Fabric Island Dressing

CUTIPLAST STERIL is a flexible fabric island dressing combining a low adherent wound pad with adhesive dressing.

## Advantages

### Absorbent, low-adherent pad

The non-woven viscose/polyester dressing pad absorbs exudate whilst the copolyester impregnated surface minimises adherence to the wound.

### Conformable

The non-woven adhesive backing conforms well to highly contoured parts of the body.

### Low allergy adhesive

The high tack polyacrylate adhesive is hypoallergenic and leaves no residue on removal.

## Indications

- Post-operative dressing, where a waterproof dressing is not required
- Unshaved areas
- Joints and other areas requiring extra flexibility
- Post minimal invasive surgery
- Post trauma
- Lacerations
- Abrasions

## Presentations

CUTIPLAST STERIL is available in a range of sizes packed in sterile peel apart packaging to enable aseptic dressing technique. CUTIPLAST is a non-sterile version and is available as a conformable fabric dressing strip in a clean, ready to use "roll" format.

### Ordering CUTIPLAST STERIL

Code	Product description (dressing; pad)	Items per unit
66001478	7.2cm x 5cm; pad: 4cm x 2.5cm	Box/100
66001473	10cm x 8cm; pad: 5.3cm x 4.5cm	Box/50
66001474	15cm x 8cm; pad: 11cm x 4cm	Box/50
66001475	20cm x 10cm; pad: 15.5cm x 5.3cm	Box/50
66001476	25cm x 10cm; pad: 20cm x 5.3cm	Box/50
66001477	30cm x 10cm; pad: 25cm x 5.3cm	Box/50

### Ordering CUTIPLAST Conformable fabric dressing strips

Code	Product description (dressing; pad)	Items per unit
66001464	4cm x 5m Roll; pad: 1.8cm x 5m (non-sterile)	Roll/1
66001466	6cm x 5m Roll; pad: 2.5cm x 5m (non-sterile)	Roll/1
66001468	8cm x 5m Roll; pad: 4cm x 5m (non-sterile)	Roll/1



Smith & Nephew Pty Ltd  
Australia  
T +61 2 9857 3999 F +61 2 9857 3900  
www.smith-nephew.com/australia  
Customer Service  
T 13 13 60 F 1800 671 000

Smith & Nephew Ltd  
New Zealand  
T +64 9 820 2840 F +64 9 820 2841  
www.smith-nephew.com/new-zealand  
Customer Service  
T 0800 657 799 (Surgical) T 0800 807 663 (Wound)

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Front cover photography copyrighted for ANZ only  
®Trademark of Smith & Nephew  
SN12744 (11/16)