



# Highly conformable retention tape

Performance with security and comfort

 **smith&nephew**  
**PRIMAFIX<sup>◊</sup> Plus**  
Conformable Retention Tape

# PRIMAFIX<sup>◊</sup> Plus, a highly conformable adhesive retention tape

PRIMAFIX Plus is a highly conformable, low allergy dressing retention tape.

PRIMAFIX Plus provides secure fixation, allowing joints to be dressed without constraining movement while providing security with patient comfort. The bandage is air and water vapour permeable allowing the skin to breathe.



## Indications

Retention of primary dressings, particularly on joints and contoured areas of the body.


## Additional application

Secure fixation of catheters and draining tubes.


## Instructions for use




**CUT**



**PULL**



**PEEL**



**PLACE**



The grid pattern facilitates measurement and cutting of the required length, ensuring easy and safe application.

## Ordering PRIMAFIX Plus

Product code	Description	Unit of Measure
36361503	PRIMAFIX Plus 2.5cm x 10m	Box/1
36361504	PRIMAFIX Plus 5cm x 10m	Box/1
36361505	PRIMAFIX Plus 10cm x 10m	Box/1
36361506	PRIMAFIX Plus 15cm x 10m	Box/1
36361509	PRIMAFIX Plus 20cm x 10m	Box/1

Smith & Nephew Pty Ltd

Australia

T 61 2 9857 3999 F 61 2 9857 3900

[www.smith-nephew.com/australia](http://www.smith-nephew.com/australia)

Customer Service

T 13 13 60 F 1800 671 000

New Zealand

T 64 9 820 2840 F 64 9 820 2841

[www.smith-nephew.com/new-zealand](http://www.smith-nephew.com/new-zealand)

Customer Service

T 0800 807 663 F 0800 263 222

<sup>◊</sup>Trademark of Smith & Nephew  
SNI2392 (05/16)