


IV3000 product range

Selection of the correct dressing for insertion sites is advised as a best practice to prevent infections¹.

IV3000 offers an extensive range of dressings designed to meet the most common IV therapy practices. Our range includes four types of dressings in various sizes²:

Ported dressings


A transparent film dressing with moisture-responsive technology³ that incorporates non-woven wings which reinforce fixation around the catheter site. The dressing has an opening designed to host the cannula's port. The dressing also includes two sterile **IMPROVED** securing strips for catheter or tube fixation and a sterile documentation label.^{4,5}



Product code	Dressing size	Quantity per carton	Recommended indication
66004011	5cm x 6cm	100	Paediatric
4006	7cm x 9cm	100	Peripheral
66004009	9cm x 12cm (oval)	50	Central
66800512	11cm x 14cm	25	PICC

Non-ported dressings


A transparent film dressing with moisture-responsive technology³ with two **IMPROVED** securing strips and documentation label.^{4,5}



Product code	Dressing size	Quantity per carton	Recommended indication
4007	6cm x 7cm	100	Peripheral
4008	10cm x 12cm	50	Central

Frame Delivery


A transparent film dressing with moisture-responsive technology.³ The frame delivery system allows user to position the dressing over the catheter with immediate visualisation of the IV site. The dressing includes two securing strips and documentation label.



Product code	Dressing size	Quantity per carton	Recommended indication
59410082	6cm x 7cm	100	Peripheral
59410882	10cm x 12cm	50	Central

Reinforced handle (orange-handles)

A transparent film dressing with moisture-responsive technology³ with edges of orange tape. These can either be torn from the edges of the dressing and used for extra security of the catheter or left in place as part of the dressing.



Product code	Dressing size	Quantity per carton	Recommended indication
4923	6cm x 8cm	100	Peripheral
4924	6cm x 8.5cm	100	Peripheral
4925	10cm x 14cm	10	Central
4973	10cm x 14cm	50	Central
4649	10cm x 20cm	50	PICC/Epidural